

Python Biax Data Reader: `biaxread`

J.R. Leeman

July 20, 2013

1 Purpose

This script was designed to read data into an easy to use format from xlook output in ASCII or binary form. Data is read, the header parsed for column names, lengths, etc., and a rec array returned that is easy to access and call from within a script. Column names are used to call the data columns, so as long as consistent naming is used the column order in the file is irrelevant.

2 Description

2.1 ASCII Files

The function ReadAscii first opens the text file with a standard open command. We know that each column in the data is written as 12 characters wide. The first line of the header is the number of records, the second is the column number, the third the column headings, the fourth the column units, and the fifth the number of records in each column. This information is parsed and stored. Numerical data is read into an array. A rec array is created with the numpy package and the labels of the column names that were parsed. Finally the array is populated with data and the final product returned to the calling script.

2.2 Binary Files

The function ReadBin opens the binary file in 'read binary' mode. The initial file information is stored as follows (in big-endian format):

Information	Format	Bytes
Name	20 characters	20
Number of Columns	int	4
Sweep	int	4
Date/Time	int	4

For each possible column there are 84 bytes of information. There are 32 possible columns. We read the information for each column and store information for columns that have data. Blank columns are identified by the first 6 characters of the name 'no_val'.

Information	Format	Bytes
Name	13 characters	13
Units	13 characters	13
Gain	int	4
Comment	50 characters	50
Number of Elements	int	4

Data is stored as the default machine format, which for modern Intel based computers is little-endian. This is the default for the module, but big-endian can be specified. Each column is written out as a sequence of doubles that are n_{elem} long. Doubles are read and stored in a numpy array. The array is then combined with datatype information collected from the headers and stored as a rec array.

To find the default machine binary format run: `python -c "import struct; print 'little' if ord(struct.pack('L', 1)[0]) else 'big'"`

A new set of improvements would be to add a cutting tool or a way to interface the script with the Pandas data handling library.

3 Cautions

A few cautions should be observed when using the `biaxread` script:

1. The empty array is shaped by row number information from the header. If a datafile is cut, but the header is left unmodified there will be many extra zero data pairs at the end of the array.

4 Usage

4.1 ASCII Files

First process the experimental data in `xlook` and output a text file with headers. To do this input *type 0 -1 1 12 pxxxx_data.txt* at the `xlook` command line or add it to the data reduction file. Be sure to correct any column or row number specifications here. Start IPython or open your own Python script. Import `biaxread` and pass the function `ReadAscii` the name of the text output file from `xlook`. The array is now returned and ready to use.

4.2 Binary Files

Process the experimental data in `xlook` and output a binary file. This is done with the 'write' command with the syntax *write filename*. In IPython or your Python script import `biaxread` and pass the function `ReadBin` the name of the binary file. The array is returned and ready to use.

Files that were written in a big-endian format by setting `dataendianness` to 'big'. The function call would look like *ReadBin(filename,dataendianness='big')*.

5 Acknowledgements

Thank you to Marco Scuderi for using the code and being patient as little issues were worked out. This code is free for all to use with proper acknowledgement. Please send any bug reports to `kd5wxb@gmail.com`.

6 Code

```
1
2 import numpy as np
3 import struct
4
5 def ReadAscii(filename):
6     """
7     Takes a filename containing the text output (with headers) from xlook and
8     reads the columns into a rec array for easy data processing and access.
9     """
10
11     try:
12         f = open(filename, 'r')
13     except:
14         print "Error Opening %s" %filename
15         return 0
16
17     col_width = 12 # Columns are 12 char wide in header
18
19     # First line of the file is the number of records
20     num_recs = f.readline()
21     num_recs = int(num_recs.strip('number of records = '))
22     print "\nNumber of records: %d" %num_recs
23
24     # Second line is column numbers, we don't care so just count them
25     num_cols = f.readline()
26     num_cols = num_cols.split('col')
27     num_cols = len(num_cols)
28     print "Number of columns: %d" %num_cols
29
30     # Third line is the column headings
31     col_headings_str = f.readline()
32     col_headings_str = col_headings_str[5:-1]
33     col_headings = []
34     for i in xrange(len(col_headings_str)/12):
35         heading = col_headings_str[12*i:12*i+12].strip()
36         col_headings.append(heading)
37
38     # Fourth line is column units
39     col_units_str = f.readline()
40     col_units_str = col_units_str[5:-1]
41     col_units = []
42     for i in xrange(len(col_units_str)/12):
43         heading = col_units_str[12*i:12*i+12].strip()
44         col_units.append(heading)
45     col_units = [x for x in col_units if x != '\n'] #Remove newlines
46
47     # Fifth line is number of records per column
48     col_recs = f.readline()
49     col_recs = col_recs.split('recs')
50     col_recs = [int(x) for x in col_recs if x != '\n']
51
52     # Show column units and headings
53     print "\n\n-----"
```

```

54     print "|%15s|%15s|%15s|" %('Name', 'Unit', 'Records')
55     print "-----"
56     for column in zip(col_headings,col_units,col_recs):
57         print "|%15s|%15s|%15s|" %(column[0],column[1],column[2])
58     print "-----"
59
60     # Read the data into a numpy recarray
61     dtype=[]
62     dtype.append(('row_num','float'))
63     for name in col_headings:
64         dtype.append((name,'float'))
65     dtype = np.dtype(dtype)
66
67     data = np.zeros([num_recs,num_cols])
68
69     i=0
70     for row in f:
71         row_data = row.split()
72         for j in xrange(num_cols):
73             data[i,j] = row_data[j]
74         i+=1
75     data_rec = np.rec.array(data,dtype=dtype)
76
77     f.close()
78
79     return data_rec
80
81 def ReadBin(filename,dataendianness='little'):
82     """
83     Takes a filename containing the binary output from xlook and
84     reads the columns into a rec array for easy data processing and access.
85
86     The data section of the file is written in the native format of the machine
87     used to produce the file. Endianness of data is little by default, but may
88     be changed to 'big' to accomodate older files or files written on power pc
89     chips.
90     """
91
92     try:
93         f = open(filename,'rb')
94     except:
95         print "Error Opening %s" %filename
96         return 0
97
98     col_headings = []
99     col_recs     = []
100    col_units     = []
101
102    # Unpack information at the top of the file about the experiment
103    name = struct.unpack('20c',f.read(20))
104    name = ''.join(str(i) for i in name)
105    name = name.split("\0")[0]
106    print "\nName: ",name
107
108    # The rest of the header information is written in big endian format

```

```

109
110 # Number of records (int)
111 num_recs = struct.unpack('>i',f.read(4))
112 num_recs = int(num_recs[0])
113 print "Number of records: %d" %num_recs
114
115 # Number of columns (int)
116 num_cols = struct.unpack('>i',f.read(4))
117 num_cols = int(num_cols[0])
118 print "Number of columns: %d" %num_cols
119
120 # Sweep (int) - No longer used
121 swp = struct.unpack('>i',f.read(4))[0]
122 print "Swp: ",swp
123
124 # Date/time(int) - No longer used
125 dtme = struct.unpack('>i',f.read(4))[0]
126 print "dtme: ",dtme
127
128 # For each possible column (32 maximum columns) unpack its header
129 # information and store it. Only store column headers of columns
130 # that contain data. Use termination at first NUL.
131 for i in range(32):
132
133     # Channel name (13 characters)
134     cname = struct.unpack('13c',f.read(13))
135     cname = ''.join(str(i) for i in cname)
136     cname = cname.split("\0")[0]
137
138     # Channel units (13 characters)
139     chunits = struct.unpack('13c',f.read(13))
140     chunits = ''.join(str(i) for i in chunits)
141     chunits = chunits.split("\0")[0]
142
143     # This field is now unused, so we just read past it (int)
144     gain = struct.unpack('>i',f.read(4))
145
146     # This field is now unused, so we just read past it (50 characters)
147     comment = struct.unpack('50c',f.read(50))
148
149     # Number of elements (int)
150     nelem = struct.unpack('>i',f.read(4))
151     nelem = int(nelem[0])
152
153     if cname[0:6] == 'no_val':
154         continue # Skip Blank Channels
155     else:
156         col_headings.append(cname)
157         col_recs.append(nelem)
158         col_units.append(chunits)
159
160
161 # Show column units and headings
162 print "\n\n-----"
163 print "|%15s|%15s|%15s|" %('Name','Unit','Records')

```

```

164     print "-----"
165     for column in zip(col_headings,col_units,col_recs):
166         print "|%15s|%15s|%15s|" %(column[0],column[1],column[2])
167     print "-----"
168
169     # Read the data into a numpy recarray
170     dtype=[]
171     for name in col_headings:
172         dtype.append((name,'double'))
173     dtype = np.dtype(dtype)
174
175     data = np.zeros([num_recs,num_cols])
176
177     for col in range(num_cols):
178         for row in range(col_recs[col]):
179             if dataendianness == 'little':
180                 data[row,col] = struct.unpack('<d',f.read(8))[0]
181             elif dataendianness == 'big':
182                 data[row,col] = struct.unpack('>d',f.read(8))[0]
183             else:
184                 print "Data endian setting invalid, please check and retry"
185                 return 0
186
187     data_rec = np.rec.array(data,dtype=dtype)
188
189     f.close()
190
191     return data_rec

```
